

JOINT ADVANCED MATERIALS & STRUCTURES
CENTER OF EXCELLENCE

Development of Dynamic Mechanical Analyzer (DMA) Calibration and Testing Procedures

2013 Technical Review

Matt Opliger

National Institute for Aviation Research,
Wichita State University

Development of Dynamic Mechanical Analyzer (DMA) Calibration and Testing Procedures

- Motivation and Key Issues
 - Determining the glass transition temperature (T_g) by DMA is critical for characterizing polymeric materials
 - Material service temperature determination
 - Quality control tool to validate curing process
 - Lack of well defined calibration and testing procedures
 - ASTM D7028 and SACMA SRM 18R refer to the equipment manufacturer's recommendations in many cases.
 - Thermocouple location
 - Temperature calibration standards
 - Leads to poor lab-to-lab and equipment-to-equipment reproducibility
 - Round robin testing of ASTM D7028 (report number D30-1004 dated December 15th, 2007) provides a poor reproducibility standard deviation of 13.36° F.
 - Significant reproducibility issues with several material qualification programs for FAA aircraft certification.
 - Accurate and reproducible T_g determination is critical to safety wherever polymer matrix composites are utilized.

Development of Dynamic Mechanical Analyzer (DMA) Calibration and Testing Procedures

- Principal Investigators & Researchers
 - John Tomblin, PhD, NIAR/WSU
 - Matthew Opliger, NIAR/WSU
 - Yeow Ng, NIAR/WSU
- Technician & Research Engineer
 - Ping Teoh
- FAA Technical Monitor
 - David Westlund
- Industry Participation
 - Joy Wu, Hexcel
 - Chuck Olson, Cytec
 - Grant Pomeroy, Intec
 - George Parker, Boeing
 - John Moylan, Delsen

Development of Dynamic Mechanical Analyzer (DMA) Calibration and Testing Procedures

- Objective
 - Reduce the lab-to-lab and equipment-to-equipment variability of Tg measurements by developing universal guidelines for temperature calibration and testing procedures for DMA equipment to enhance pre-existing testing standards.
 - Improve industry applications and safety
 - Material service temperature definition
 - Quality control

Development of Dynamic Mechanical Analyzer (DMA) Calibration and Testing Procedures


- Approach
 - Research areas that potentially influence Tg determination and reproducibility.
 - Thermocouple location
 - Specimen dimensions
 - Temperature calibration
 - Calibration materials
 - Application of calibration materials
 - Develop universal guidelines for temperature calibration and testing procedures for DMA equipment.
 - Evaluate developed guidelines for lab-to-lab and equipment-to-equipment reproducibility through round robin testing.
 - Further evaluate developed guidelines by comparing to DSC and TMA test results.

Thermocouple Location


- Thermocouple (TC) Location Evaluation
 - ASTM D7028 section 11.9 remarks, “*Follow the manufacturer’s manual or recommendations to position the thermocouple. Typically the thermocouple should be as close to the sample as possible.*”
 - Manufacturers recommend that the thermocouple proximity should be close to the sample, as well, but do not provide any specifics.
 - The fixture acts as a heat sink, affecting heat transfer to the sample.
 - Evaluate TC location, relative to location to fixture
 - Evaluate TC location, relative to location to sample

Thermocouple Location – Test Configurations

TC Location	Test Method	Specimens per Material (Number of Materials)		Material ID
		50mm 3-pt Bend Fixture	35mm Dual Cantilever Fixture	
1	ASTM D7028 & Figure 1	1 x (4)	1 x (4)	180F-Tg-EPXY 350F-Tg-EPXY 500F-Tg-BMI 600F-Tg-PI
2		1 x (4)	1 x (4)	180F-Tg-EPXY 350F-Tg-EPXY 500F-Tg-BMI 600F-Tg-PI
3		1 x (4)	1 x (4)	180F-Tg-EPXY 350F-Tg-EPXY 500F-Tg-BMI 600F-Tg-PI
4		1 x (4)	1 x (4)	180F-Tg-EPXY 350F-Tg-EPXY 500F-Tg-BMI 600F-Tg-PI
5		1 x (4)	1 x (4)	180F-Tg-EPXY 350F-Tg-EPXY 500F-Tg-BMI 600F-Tg-PI
6		1 x (4)	1 x (4)	180F-Tg-EPXY 350F-Tg-EPXY 500F-Tg-BMI 600F-Tg-PI




TC1




TC3


Top



Side




TC6



Thermocouple Location – Test Results


35mm Dual Cantilever Fixture



Approximate Location of Furnace Wall

Thermocouple Location – Test Results

50mm 3-Point Bend Fixture



Approximate Location of Furnace Wall

Thermocouple Location – Conclusion

- Thermocouple (TC) Location Evaluation
 - Significant differences between TC locations
 - The fixture acts as a heat sink, affecting heat transfer to the sample. This is evident by low T_g values observed at TC1 and TC2.
 - The proximity of the test specimen to the furnace becomes more critical at higher temperatures
- Recommend TC3 Location
 - Locations TC2, TC4, and TC6 can result in a damaged TC if a highly deformable material is tested.
 - Location TC3 is least affected by the fixture heat sink and furnace proximity.
 - T_g values obtained at TC3 were in the top half of all TC locations


Specimen Dimensions

- Specimen Dimensions Evaluation
 - Determine if specimen dimensions affect T_g
 - Evaluate changes in thickness
 - Evaluate changes in width
 - Evaluate changes in length
 - Possible effects include
 - Thermal lag
 - Specimen stiffness

Specimen Dimensions – Test Configurations

Property	Test Method	Support Dimensions			Specimens per Material (Number of Materials)		Material ID
		Thickness (in.)	Width (in.)	Length (in.)	50mm 3-pt Bend Fixture	35mm Dual Cantilever Fixture	
Glass Transition Temperature (T _g), Dry, by DMA	ASTM D7028 & Figure 2	0.08	0.5	2.5	1 x (1)	1 x (1)	350F-Tg-EPXY
		0.12	0.5	2.5	1 x (1)	1 x (1)	350F-Tg-EPXY
		0.04	0.5	2.5	1 x (1)	1 x (1)	350F-Tg-EPXY
		0.08	0.25	2.5	1 x (1)	1 x (1)	350F-Tg-EPXY
		0.08	0.5	2.25	1 x (1)	1 x (1)	350F-Tg-EPXY

Figure 2



Thickness X Width X Length
0.08" X 0.5" X 2.5"
a. Baseline Specimen

0.12" X 0.5" X 2.5"
b. Variable Thickness Specimen

0.04" X 0.5" X 2.5"
c. Variable Thickness Specimen


0.04" X 0.5" X 2.5"
d. Variable Width Specimen

0.08" X 0.5" X 3.0"
e. Variable Length Specimen

Specimen Dimensions – Test Results

Property	Test Method	Support Dimensions			T _g [°F]		Material ID
		Thickness (in.)	Width (in.)	Length (in.)	50mm 3-pt Bend Fixture	35mm Dual Cantilever Fixture	
Glass Transition Temperature (T _g), Dry, by DMA	ASTM D7028 & Figure 2	0.08	0.5	2.5	371	374	350F-Tg-EPXY
		0.12	0.5	2.5	365	375	350F-Tg-EPXY
		0.04	0.5	2.5	385	384	350F-Tg-EPXY
		0.08	0.25	2.5	381	385	350F-Tg-EPXY
		0.08	0.5	2.25	372	373	350F-Tg-EPXY

Figure 2



Thickness X Width X Length
0.08" X 0.5" X 2.5"
a. Baseline Specimen

0.12" X 0.5" X 2.5"
b. Variable Thickness Specimen

0.04" X 0.5" X 2.5"
c. Variable Thickness Specimen

0.04" X 0.5" X 2.5"
d. Variable Width Specimen




0.08" X 0.5" X 3.0"
e. Variable Length Specimen

Specimen Dimensions – Conclusion




Thickness dependent material behavior

0.12" (Thick) ———
 0.04" (Thin) - - - -

50mm 3-Pt Bend Fixture



35mm Dual Cantilever Fixture




- Temperature Calibration Evaluation
 - ASTM D7028, section 9.1 states, “*The DMA equipment shall be calibrated in accordance with Test Method E 1867 for temperature signals and in accordance with the equipment manufacturer’s recommendation for the storage modulus.*”
 - ASTM E1867 lists NIST traceable temperature calibration standards
 - Indium, $T_m = 156.6^\circ \text{C}$
 - Tin, $T_m = 231.9^\circ \text{C}$
 - Lead, $T_m = 327.5^\circ \text{C}$
 - Zinc, $T_m = 419.5^\circ \text{C}$
 - ASTM E 1867 does not address the need for a rigid span support to hold the calibration standard.
 - Evaluate calibration support beam dimensions/configuration
 - Evaluate calibration support beam materials

Temperature Calibration – Support Dimensions & Configurations

Property	Test Method	Support Dimensions			Specimens per Material (Number of Materials)		Material ID
		Thickness (in.)	Width (in.)	Length (in.)	50mm 3-pt Bend Fixture	35mm Dual Cantilever Fixture	
Indium Melting Temperature	ASTM E1867 & Figure 3	0.04	0.5	Bottom, 2.5 Top, 2.5	1 x (1)	1 x (1)	500F-Tg-BMI
		0.02	0.5	Bottom, 2.5 Top, 2.5	1 x (1)	1 x (1)	500F-Tg-BMI
		0.06	0.5	Bottom, 2.5 Top, 2.5	1 x (1)	1 x (1)	500F-Tg-BMI
		0.04	0.25	Bottom, 2.5 Top, 2.5	1 x (1)	1 x (1)	500F-Tg-BMI
		0.04	0.5	Bottom, 2.5 Top, 1.5	1 x (1)	1 x (1)	500F-Tg-BMI
Tin Melting Temperature		0.04	0.5	Bottom, 2.5 Top, 2.5	1 x (1)	1 x (1)	500F-Tg-BMI
		0.02	0.5	Bottom, 2.5 Top, 2.5	1 x (1)	1 x (1)	500F-Tg-BMI
		0.06	0.5	Bottom, 2.5 Top, 2.5	1 x (1)	1 x (1)	500F-Tg-BMI
		0.04	0.25	Bottom, 2.5 Top, 2.5	1 x (1)	1 x (1)	500F-Tg-BMI
		0.04	0.5	Bottom, 2.5 Top, 1.5	1 x (1)	1 x (1)	500F-Tg-BMI

Figure 3




Bottom Thickness, Top Thickness X Bottom Width, Top Width X Bottom Length, Top Length
 0.04", 0.04" X 0.5", 0.5" X 2.5", 2.5" 0.02", 0.02" X 0.5", 0.5" X 2.5", 2.5" 0.06", 0.06" X 0.5", 0.5" X 2.5", 2.5" 0.04", 0.04" X 0.25", 0.25" X 2.5", 2.5" 0.04", 0.04" X 0.5", 0.5" X 2.5", 1.5"
 a. Baseline Specimen b. Variable Thickness Specimen c. Variable Thickness Specimen d. Variable Width Specimen e. Variable Length Specimen



Temperature Calibration – Support Dimensions & Configurations, Test Results

Property	Test Method	Support Dimensions			T _m [°F]		Material ID
		Thickness (in.)	Width (in.)	Length (in.)	50mm 3-pt Bend Fixture	35mm Dual Cantilever Fixture	
Indium Melting Temperature (313.88°F)	ASTM E1867 & Figure 3	0.04	0.5	Bottom, 2.5 Top, 2.5	323	335	500F-Tg-BMI
		0.02	0.5	Bottom, 2.5 Top, 2.5	321	333	500F-Tg-BMI
		0.06	0.5	Bottom, 2.5 Top, 2.5	323	334	500F-Tg-BMI
		0.04	0.25	Bottom, 2.5 Top, 2.5	326	344	500F-Tg-BMI
		0.04	0.5	Bottom, 2.5 Top, 1.5	321	334	500F-Tg-BMI
Tin Melting Temperature (449.38°F)	ASTM E1867 & Figure 3	0.04	0.5	Bottom, 2.5 Top, 2.5	458	471	500F-Tg-BMI
		0.02	0.5	Bottom, 2.5 Top, 2.5	455	467	500F-Tg-BMI
		0.06	0.5	Bottom, 2.5 Top, 2.5	458	467	500F-Tg-BMI
		0.04	0.25	Bottom, 2.5 Top, 2.5	458	472	500F-Tg-BMI
		0.04	0.5	Bottom, 2.5 Top, 1.5	456	470	500F-Tg-BMI

Figure 3



Bottom Thickness, Top Thickness X Bottom Width, Top Width X Bottom Length, Top Length

0.04", 0.04" X 0.5", 0.5" X 2.5", 2.5"

a. Baseline Specimen

0.02", 0.02" X 0.5", 0.5" X 2.5", 2.5"

b. Variable Thickness Specimen

0.06", 0.06" X 0.5", 0.5" X 2.5", 2.5"

c. Variable Thickness Specimen

0.04", 0.04" X 0.25", 0.25" X 2.5", 2.5"

d. Variable Width Specimen

0.04", 0.04" X 0.5", 0.5" X 2.5", 1.5"

e. Variable Length Specimen

Temperature Calibration – Support Material

Polymer Matrix Material for Comparison to Evaluation of Support Materials

Property	Test Method	Specimens per Material (Number of Materials)		Material ID
		50mm 3-pt Bend Fixture	35mm Dual Cantilever Fixture	
Indium Melting Temperature	ASTM E1867	3 x (3)	3 x (3)	350F-Tg-EPXY 500F-Tg-BMI 600F-Tg-PI
Tin Melting Temperature		3 x (2)	3 x (2)	500F-Tg-BMI 600F-Tg-PI
Specific Heat Capacity, Thermal Conductivity, and Thermal Diffusivity	ASTM E1461	1 x (4)		180F-Tg-EPXY 350F-Tg-EPXY 500F-Tg-BMI 600F-Tg-PI

Support Material Evaluation

Property	Test Method	Specimens per Material (Number of Materials)		Material ID
		50mm 3-pt Bend Fixture	35mm Dual Cantilever Fixture	
Indium Melting Temperature	ASTM E1867	3 x (4)	3 x (4)	Macor Alumina Silicate Silica Glass Steel
Tin Melting Temperature		3 x (4)	3 x (4)	Macor Alumina Silicate Silica Glass Steel
Specific Heat Capacity, Thermal Conductivity, and Thermal Diffusivity	ASTM E1461	1 x (4)		Macor Alumina Silicate Silica Glass Steel

Temperature Calibration – Support Material, Test Results

Indium Tm = 313.88°F, Tin Tm = 449.38°F

50mm 3-pt Bend Fixture		
Material ID	Indium Tmavg [°F]	Tin Tmavg [°F]
350F-Tg-EPXY	319.77	-
500F-Tg-BMI	320.06	454.42
600F-Tg-PI	317.72	453.16
Polymer Matrix Avg	319.18	453.79
Alumina Silicate	322.84	455.74
Macor	316.74	453.01
Silica Glass	317.61	453.72
Steel	330.83	461.62

35mm Dual Cantilever Fixture		
Material ID	Indium Tmavg [°F]	Tin Tmavg [°F]
350F-Tg-EPXY	330.27	-
500F-Tg-BMI	330.41	463.80
600F-Tg-PI	327.84	462.52
Polymer Matrix Avg	329.51	463.16
Alumina Silicate	330.34	463.79
Macor	328.50	460.49
Silica Glass	328.16	460.80
Steel	336.60	467.14

Temperature Calibration – Support Material, Test Results

Material	Thermal Diffusivity [mm^2/s]		
	α at T = 77°F	α at T = 300°F	α at T = 450°F
180F-TG-EPOXY	0.255		
350F-TG-EPOXY	0.407	0.331	
500F-TG-BMI	0.404	0.341	0.310
600F-TG-PI	0.418	0.348	0.319

Material	Description	Thermal Diffusivity [mm^2/s]		
		α at T = 77°F	α at T = 300°F	α at T = 450°F
Macor	Glass-Ceramic	0.844	0.727	0.684
Lavastone	Alumina-Silicate	2.00	1.39	1.17
Soda Lime Glass	Silicate Glass	0.505	0.477	0.461

Temperature Calibration – Conclusion

- Calibration Support Beam Geometry
 - No significant differences seen with variable dimensions
 - Recommend using a similar size calibration support beam as the intended test specimen
- Calibration Support Beam Material
 - Macor, Lavastone and soda lime glass were all found to be good materials.
 - Comparable melting temperatures to melting temperatures found using polymer matrix composite calibration support beams
 - Similar thermal diffusivities to polymer matrix composites
 - Recommend using Macor
 - Relatively easy to obtain
 - Easily machineable
 - Consistent formulation

Development of Dynamic Mechanical Analyzer (DMA) Calibration and Testing Procedures

- Development of Calibration and Testing Procedures
 - Utilized data from the DMA evaluations
 - TC location
 - Specimen dimensions
 - Span support dimensions/configuration
 - Span support material
 - Distributed new calibration and testing procedures and materials to labs for round robin testing for evaluation of procedures.
 - Low-to-high T_g materials are being evaluated
 - Specimens are being tested after dry and wet environmental conditioning

Inter-Laboratory Testing – Test Matrix

Property	Test Method	Specimens per Material (Number of Materials)												Material ID
		Laboratory, 50mm 3-pt Bend Fixture						Laboratory, 35mm Dual Cantilever Fixture						
		Lab 1	Lab 2	Lab 3	Lab 4	Lab 5	Lab 6	Lab 1	Lab 2	Lab 3	Lab 4	Lab 5	Lab 6	
DMA Temperature Calibration	ASTM E1867	1 x (2)	1 x (2)	1 x (2)	1 x (2)	1 x (2)	1 x (2)	1 x (2)	1 x (2)	1 x (2)	1 x (2)	1 x (2)	1 x (2)	Indium Zinc
DMA Temperature Verification		1 x (2)	1 x (2)	1 x (2)	1 x (2)	1 x (2)	1 x (2)	1 x (2)	1 x (2)	1 x (2)	1 x (2)	1 x (2)	1 x (2)	Indium Zinc
Glass Transition Temperature (T _g), Dry, by DMA	ASTM D7028	3 x (5)	3 x (5)	3 x (5)	3 x (5)	3 x (5)	3 x (5)	3 x (5)	3 x (5)	3 x (5)	3 x (5)	3 x (5)	3 x (5)	180F-Tg-EPXY 350F-Tg-EPXY 500F-Tg-BMI 600F-Tg-PI Polycarbonate
Glass Transition Temperature (T _g), Wet, by DMA		3 x (4)	3 x (4)	3 x (4)	3 x (4)	3 x (4)	3 x (4)	3 x (4)	3 x (4)	3 x (4)	3 x (4)	3 x (4)	3 x (4)	180F-Tg-EPXY 350F-Tg-EPXY 500F-Tg-BMI 600F-Tg-PI

Inter-Laboratory Testing – Test Results

- Round Robin Test Results
 - Tests have been completed from three of the six labs
 - The standard deviation of the dry tested samples is much improved from the ASTM D 7028 round robin.
 - The standard deviation of the wet tested samples is worse than that obtained from the ASTM D 7028 round robin.
 - Two of the labs required relatively thin (0.06”) samples and the third and remaining labs are using a relatively thick (0.12”) sample.
 - There is reason to believe that the thin samples lost more moisture than the thick samples before the Tg was reached.
 - All round robin results will be reported in the FAA report.

	Standard Deviation of the Glass Transition Temperature					
	[°F]					
	Dry T _{gE'}	Dry T _{gE''}	Dry T _{g_{stand}}	Wet T _{gE'}	Wet T _{gE''}	Wet T _{g_{stand}}
ASTM D 7028 Round Robin (2007)	18.628	18.753	16.259	13.254	12.544	13.281
Current Round Robin w/ Calibration Procedures	6.640	6.568	5.325	15.476	17.129	21.750

Development of Dynamic Mechanical Analyzer (DMA) Calibration and Testing Procedures


- Viability of New Calibration and Testing Procedures
 - Data from round robin testing will be analyzed to determine adequacy of the new procedures.
 - Additional thermal analysis techniques were used on a small scale to further characterize the capabilities of these procedures.
 - Differential Scanning Calorimetry (DSC) by ASTM D3418
 - Thermo-Mechanical Analysis (TMA) by expansion method ASTM E1545 (procedure A)
 - Thermo-Mechanical Analysis (TMA) by flexure method ASTM E1545 (modified procedure B)

Comparison to DSC & TMA Techniques



Property	Test Method	Specimens per Material (Number of Materials)	Material ID
Glass Transition Temperature (T _g), Dry, by DSC	ASTM D3418	3 x (5)	180F-Tg-EPXY 350F-Tg-EPXY 500F-Tg-BMI 600F-Tg-PI Polycarbonate
Glass Transition Temperature (T _g), Dry, by TMA Expansion	ASTM E1545, Procedure A	3 x (5)	180F-Tg-EPXY 350F-Tg-EPXY 500F-Tg-BMI 600F-Tg-PI Polycarbonate
Glass Transition Temperature (T _g), Dry, by TMA Flex	ASTM E1545, Modified Procedure B	3 x (5)	180F-Tg-EPXY 350F-Tg-EPXY 500F-Tg-BMI 600F-Tg-PI Polycarbonate

Comparison to DSC & TMA Techniques – Test Results



Material	Specimen Number	Glass Transition Temperature [°F]						
		DMA by ASTM D 7028 (NIAR Calibration Procedures)			TMA Expansion by ASTM E 1545 (Procedure A)	TMA Flexure by ASTM E 1545 (Procedure B)	DSC by ASTM D 3418	
		Dry T _{gE'}	Dry T _{gE''}	Dry T _{gTanδ}	Dry T _{gExp}	Dry T _{gFlex}	Dry T _{gEig}	Dry T _{gmg}
180F-TG-EPOXY	Mean	142	206	214	149	176	154	165
350F-TG-EPOXY	Mean	363	375	388	364	361	363	375
500F-TG-BMI	Mean	519	548	559	520	503	506	527
600F-TG-PI	Mean	614	641	657	609	609	601	618
Polycarbonate	Mean	297	298	313	301	297	291	301



Comparison to DSC & TMA Techniques – Test Results




Comparison to DSC & TMA Techniques – Test Results



Looking Forward


- Benefit to Aviation
 - Improvement in industry applications and safety
 - Material service temperature definition
 - Quality control
- Future needs
 - Further evaluate differences observed due to dimensional variation of specimens
 - Look at equations used to determine the elastic and viscous modulus
 - Look at potential thermal lag issues
 - Study methods to improve wet Tg reliability

Looking Forward – Moisture Loss During Heat Up



End of Presentation.

Thank you.



JOINT ADVANCED MATERIALS & STRUCTURES
CENTER OF EXCELLENCE

