



Molecular Diagnostics Laboratory

WICHITA STATE UNIVERSITY

In Person Specimen Drop-off/Test Kit Pick-up Instructions

WHO:

Your organization must ensure that whomever you assign to drop-off specimens and pick-up supplies for your organization follows all procedures below and has read and understands the [Category B Infectious Substance Transportation](#) document and process covered in training.

WHAT:

Specimens must be sealed individually in the provided biohazard bags and placed in a rigid container with a biohazard sticker for transport.

WHEN:

Specimen Drop-off/Test Kit Pick-up Hours

Monday – Thursday, 7:00am – 6:00pm

Friday, 7:00am – 2:00pm

WHERE:

4174 S Oliver, Building 174H, Door #2

Southeast of MacArthur and Oliver in the Spirit AeroSystems complex.

Press the video doorbell button. When you hear the click, the door is unlocked.

WHY:

Please make sure all orders are submitted & released through TestDirectly **BEFORE** dropping off the physical specimens at the lab. If a specimen is received without a released complete and accurate order, it will be rejected/discarded requiring that individual to re-test.

12/5/2023